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No Risk of Cell Damage Due to Flow Influences When Using the Flowmax® 242i Generation 3

This conclusion was reached by ReGEN Biomedical on behalf of MIB GmbH / Flowmeister. The Flowmax® 242i is a highly advanced flow meter that precisely measures the volumetric flow rate of liquids. Based on the ultrasonic technology, it is able to measure non-invasively both conductive and non-conductive liquids. All parts in contact with the media are made of PELD. The sensor has USP Class VI certification and is resistant to gamma radiation, making it excellent for single-use applications not only in the pharmaceutical industry. Due to the integrated microchip, the calibration data of the sensor are stored in the main body, making the Flowmax® 242i a 'plug and play' solution.

The technical specifications of the Flowmax® 242i allow the measurement of a wide range of liquids, including DI water, bioprocess fluids, cosmetics, liquid foods, as well as aggressive media or blood and blood-like substances. The task assigned to ReGEN Biomedical was to investigate whether the use of the Flowmax® 242i could lead to cell damages.

Simulations using COMSOL showed no turbulence; instead, a laminar flow was observed through the sensor. This is particularly important because turbulent flows can potentially have harmful effects on cells. In simulated tests with water and blood-like fluids, the possible shear forces remained below the established threshold of 0.01 Pa – a benchmark for stress on endothelial cells in the lymphatic system.

The theoretical results were experimentally validated with fibroblasts, where no reduction in cell viability was observed. This confirmation is crucial as it suggests that the Flowmax® sensor does not have harmful effects on human cells under the tested conditions. In particular, the laminarity of the flow and low shear forces indicate that the sensor can be safely used in applications requiring interaction with living cells.

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